



BOECO Pipette Fillers Balls, natural rubber, red, with 3 valves

Cat. No.	Description
BOE 100	Model STANDARD, short cone For pipettes till 25 ml
BOE 101	Model UNIVERSAL, long cone for all pipettes



BOECO Pipette Pumps, PP, can be operated with one hand. Glass and plastic pipettes fit easily. Precise filling or emptying is done by rotating their actuator-wheels forward or backward. For rapid emptying simply depress plunger.

Cat. No.	Description
BOE 120	0 - 2 ml, blue
BOE 121	0 - 10 ml, green
BOE 122	0 - 25 ml, red



BOECO Macro Pipetting Aid offers safer, more convenient filling and dispensing, with glass and plastic pipettes. Vacuum is applied by an easypull slide lever, eliminating the hand stress caused by conventional "squeeze bulb" apparatus. Simply pull on the slide whenever vacuum is needed. Onehanded control of aspiration, dispensing and vacuum saves time and leaves the other hand free to manipulate samples. The valve is in-line filtersystem, combined with the lightweight handle, provides outstanding operation, control, and safety for manual pipetting procedures.

- * A single rocker-switch accurately controls aspiration and delivery
- * Disposable, in-line filter prevents contamination
- * Holds 0.1 to 100ml pipets
- * Filter and pipet connector are autoclavable for decontamination

Cat. No.	Description
BOE 150	BOECO Macro Pipetting Aid, complete
Accessories	
BOE 150-1	Filter disc, 3µm, pack of 10
BOE 150-2	Asperating cone
BOE 150-3	Silicon adapter



BOECO Electronic Macro Pipetting Aid is a cordless rechargeable pipetting aid which can be used with all plastic or glass pipets between 1-100 ml, pasteur pipettes and optional also with 5 ml and 10 ml Pipettor tips. The linear speed control allows the aspiration and dispensing to be adjusted easily and precisely to suit either large, or small volume pipets. The speed can also be fine-tuned by applying varying finger pressure to the operation buttons.

- * The attached stand on the unit supports unit and pipet when not in use.
- * The unit is equipped with a low battery indicator

Cat. No.	Description
BOE 9500000	BOECO Electronic Macro Pipetting Aid, complete with AC-adapter 220V, and 0,45µm hydrophobic filters
BOE 9500001	Idem with AC-adapter 110V, U.S. Plug
Accessories	
BOE 712912	Replacement filter, nonsterile (5 pcs.)
BOE 712728	Nose cone and silicon adapter, 5 ml
BOE 712729	Nose cone and silicon adapter, 10 ml
BOE 985000	5 ml Pipette tip (bag of 300 pcs.)
BOE 980320	10 ml Pipette tip (100 pcs. box)